

# Westernblot

## Solutions and reagents:

### lysis buffers:

#### *RIPA buffer (radioimmunoprecipitation assay buffer):*

50 mM Tris-HCl, pH 8.0  
150 mM NaCl  
1% NP-40 or 0.1% Triton X-100  
0.5% sodium deoxycholate  
0.1% SDS  
protease inhibitors

#### *NP-40 buffer:*

50 mM Tris-HCl, pH 8.0  
150 mM NaCl  
1.0% NP-40 (or 0.1% Triton X-100)  
protease inhibitors

### 3 x Laemmli buffer/ sample buffer:

150 mM Tris-HCl, pH 6.8  
300 mM DTT  
6% SDS  
0.3% bromophenol blue  
30% glycerol

### 10x PBS:

1.37 M NaCl  
0.027 M KCl  
0.1 M Na<sub>2</sub>HPO<sub>4</sub>  
0.018 M KH<sub>2</sub>PO<sub>4</sub>

- Resolve in 800 ml ddH<sub>2</sub>O.
- Adjust pH to 7.4 using HCl.
- Fill it up to 1l.
- Autoclave it.
- Store at room temperature.
- Dilute 1:10 before use.

### blocking buffer:

3-5% non-fat dry milk or BSA in PBST (PBS + 0.1% Tween 20)

## **Procedure**

### **Sample preparation**

#### **Preparation of lysate from cell culture:**

- Place the cell culture dish on ice and wash the cells with ice-cold PBS.
- Aspirate the PBS and add ice-cold lysis buffer (1 ml per 10 cm dish).
- Scrape adherent cells off the dish using a cold plastic cell scraper and gently transfer the cell suspension into a precooled microcentrifuge tube.
- If required, cells can be harvested by trypsinization and washed with PBS prior to resuspension in lysis buffer.
- Incubate at 4°C for 30 min with constant agitation.
- Centrifuge at 16,000 x g for 20 min at 4°C.
- Transfer the supernatant to a fresh tube on ice, and discard the pellet.
- Remove a small volume (10-20 µl) of lysate for analysis by a protein assay. Determine the protein concentration for each cell lysate.
- If necessary, aliquot the protein samples for long-term storage at -20°C. Repeated freeze and thaw cycles cause protein degradation and should be avoided.
- Add ½ volume of 3x Laemmli sample buffer.
- Boil each cell lysate in sample buffer at 95°C for 5 min.
- Centrifuge at 16,000 x g in a microcentrifuge for 1 min.

### **Protein separation by SDS-PAGE**

Polyacrylamid percentage of SDS-gel for best resolution of proteins based on their molecular weight:

Protein size	Gel percentage
4-40 kDa	20%
12-45 kDa	15%
10-70 kDa	12.5%
15-100 kDa	10%
25-200 kDa	8%

- Load equal amounts of protein into the wells of the SDS-PAGE (10-50 µg/lane protein of cell lysate or 10-100 ng/lane purified protein). Add molecular weight marker in one of the lanes.
- Run the gel according to manufacturer's instructions (e.g. 1-2 h at 200 V).

### **Protein transfer from gel to membrane**

Use either nitrocellulose or PVDF membrane. Activate PVDF with methanol for 1 min and rinse with transfer buffer before preparing the stack.

Follow the manufacturer's instructions for blotting.

### **Antibody incubation**

- After transfer briefly rinse the membrane in distilled water or PBST.
- Block the membrane 1 h at RT with blocking buffer.
- Incubate the membrane with appropriate dilutions of primary antibody in blocking buffer ON at 4°C or 1-2 h at RT.
- Wash the membrane with PBST for 15 min/ 3x 5 min.
- Incubate the membrane with recommended dilution of conjugated secondary antibody in blocking buffer for 1 h at RT.
- Wash the membrane with PBST for 15 min/ 3x 5 min.
- For signal development follow the kit manufactures instructions of the detection kit used.