

PROGEN
passion for research

RELAXED
IHC
HANDLING



prohisto
IHC staining kits

EASY HANDLING KITS FOR BETTER IMMUNOHISTOCHEMISTRY



EASY HANDLING KITS

prohisto-HiSec + prohisto-HiVis


anti-mouse HRP
anti-guinea pig HRP
anti-rabbit HRP

DAB duo

Enjoy premium performance and care-free handling of PROGEN's IHC staining kits that supplement our comprehensive portfolio of primary Abs!

All HiSec antibodies also available as sets with prohisto-HiVis DAB duo.

EXTERNALLY VALIDATED FOR IHC



A) IHC of human skin using rabbit anti-Histidine decarboxylase polyclonal serum and **prohisto-HiSec anti-rabbit HRP**.

B) IHC of human colon carcinoma using guinea pig anti-Keratin K8/K18 antibody and **prohisto-HiSec anti-guinea pig HRP**, visualized with **prohisto-HiVis DAB duo**.

C) IHC of rat colon using mouse anti-Keratin K20 antibody and **prohisto-HiSec anti-mouse HRP**, visualized with **prohisto-HiVis DAB duo**.

All images courtesy of J. Heß, University Hospital Heidelberg

SHOP ONLINE at www.progen.com

& check out our IHC staining solutions!

